



INSTRUMENTS

DENTASEPT®

Special Rotative

INDICATIONS

Ready to use detergent disinfectant solution for the treatment of burs, rotary and canal instruments, prior to sterilisation.

ADVANTAGES

- ✓ Ready to use solution
- ✓ Compatible with a wide range of materials, nickel, titanium, tungsten carbide, aluminium, etc
- ✓ Contact time: 15 minutes
- ✓ Excellent microbiological properties - effective on all enveloped viruses¹
- ✓ Can be used in an Ultrasonic cleaner
- ✓ Medical Device Class IIb (Directive 93/42/EEC modified)

¹ The Robert Koch Institute recommends the use of disinfectants claiming limited virucidal activity to interrupt the infection-chain. According to the RKI, limited virucidal activity (effective against enveloped viruses) corresponds to an activity against Vaccinia virus and BVDV.

MICROBIOLOGICAL PROPERTIES AT 20°C

Bactericide: EN 1040, EN 13727+A1, EN 14561 - 5 minutes.
Yeasticide: EN 1275, EN 13624, EN 14562 - 15 minutes.
Tuberculocide: EN 14348, EN 14563 - 60 minutes.
Active in 5 minutes against PRV (surrogate of HBV), Herpes virus, BVDV (surrogate of HCV) and Vaccinia virus.

QUALITATIVE COMPOSITION

N,N-didecyl-N,N-dimethylammonium carbonate, non-ionic surfactants, sequestering agent, fragrance, colouring agent, excipients.

PACKAGING

2 litre container	2480.062
5 litre container	2480.047



INSTRUMENTS

DENTASEPT® Special Rotative PROTOCOL FOR USE

Ready to use detergent disinfectant solution for the treatment of burs, rotary and canal instruments, prior to sterilisation.

BEFORE USE

Wear appropriate gloves.

HOW TO USE

- 1 Pour Dentasept® Special Rotative in to a soaking tank. Do not dilute.
- 2 Totally immerse the instruments in to the solution, cover and leave to soak. Respect a contact time of 15 minutes.
- 3 Remove instruments from the solution. Rinse under running water.
- 4 Dry instruments properly and prepare for sterilisation.

Change the solution at least once daily.

PRECAUTIONS

Dangerous - respect the precautions for use. Medical Device Class IIb, CE 0459. Before use, read the label and information concerning the product.

Store between +5°C and +35°C.

