

Technical Datasheet

Disinfection wipes Hygowipe | PET

Art.no.: 31553

HYGO  CLEAN

Product image



Product specifications

- material: 100% PET, cloth quality with textile character
- individual impregnation (surface disinfection or cleaning substance)
- very good absorption behavior
- economical in use, without scattering loss

Specification

Colour	white
Length	36cm
Width	20cm
Grammage	40gr-m ²
tearoffs	90

Product properties

Properties	on roll
material	PET

Technical Datasheet

Disinfection wipes Hygowipe | PET

Art.no.: 31553



Packing

Type of inner packaging	PE
Type of outer packaging	cardboard
Quantity inner packing	1 roll
Quantity outer packing	6 roll
Quantity CCG1 pallet	288 roll
Quantity full pallet	648 roll
LxWxH outer packing in cm	32,0 x 22,0 x 22,5
Volume outer packing in m ³	0.01584
Weight outer packing in gr.	1.830

Additional information

Costums tariff number	6307103000
Country of origin	CN
EAN code outer packing	4015544915532
EAN code inner packing	4015544315684

You need more detailed information? We will be happy to advise you!

The technical data are average production values and may deviate slightly in individual cases.
Subject to change without notice, no warranty for misprints or errors.