

“STERIDIAMOND” STERILIZATION FLAT REELS

- Flat rolls imprints are performed outside the packing area, along the welding line, between paper and film, to prevent ink from migrating to the product.
- Transparent green film grants the prompt identification of the content.
- The film is reinforced in order to avoid any tear upon roll opening.
- High quality “medical grade” paper with 60 or 70 gsm.
- Three-line sealing with strong strength and effective spacing grants a perfect adhesivity during sterilization process and an excellent fibre-free and tear-free peeling during opening.
- Steam and EO process indicators.
- Formaldehyde process indicator is available upon request.
- ISO 11140 certified inks are water based and non-toxic.
- All materials comply with international standards.
- Opening direction is clearly marked.
- Validated production process.
- Each roll is wrapped in a protective film and placed in the packing carton to ensure protection against dust even after opening of the box.
- Suitable to be used with all types of sealing machines.

STERIDIAMOND FLAT REELS 100 M



STERIDIAMOND FLAT REELS 200 M



SRP - STERIDIAMOND FLAT REELS 100 M

REF.	SIZE	MIN LOT	PACK
SRP0510	55 mm X 100 m	8 pcs	8 pcs
SRP7510	75 mm X 100 m	8 pcs	8 pcs
SRP1010	100 mm X 100 m	4 pcs	4 pcs
SRP1510	150 mm X 100 m	4 pcs	4 pcs
SRP2010	200 mm X 100 m	2 pcs	2 pcs
SRP2510	250 mm X 100 m	2 pcs	2 pcs
SRP3010	300 mm X 100 m	2 pcs	2 pcs
SRP3510	350 mm X 100 m	1 pcs	1 pcs
SRP4010	400 mm X 100 m	1 pcs	1 pcs
SRP5010	500 mm X 100 m	1 pcs	1 pcs

SRP - STERIDIAMOND FLAT REELS 200 M

REF.	SIZE	MIN LOT	PACK
SRP0520	55 mm X 200 m	8 pcs	8 pcs
SRP7520	75 mm X 200 m	8 pcs	8 pcs
SRP1020	100 mm X 200 m	4 pcs	4 pcs
SRP1520	150 mm X 200 m	4 pcs	4 pcs
SRP2020	200 mm X 200 m	2 pcs	2 pcs
SRP2520	250 mm X 200 m	2 pcs	2 pcs
SRP3020	300 mm X 200 m	2 pcs	2 pcs
SRP3520	350 mm X 200 m	1 pcs	1 pcs
SRP4020	400 mm X 200 m	1 pcs	1 pcs




“STERIDIAMOND” STERILIZATION GUSSETTED REELS

- Manufactured with a gusset that allows bulky objects to be placed inside
- Gussetted rolls imprints are performed outside the packing area, along the welding line, between paper and film, to prevent ink from migrating to the product.
- Transparent green film grants the prompt identification of the content.
- The film is reinforced in order to avoid any tear upon roll opening.
- High quality “medical grade” paper with 60 or 70 gsm.
- Three-line sealing with strong strength and effective spacing grants a perfect adhesivity during sterilization process and an excellent fibre-free and tear-free peeling during opening.
- Steam and EO process indicators.
- Formaldehyde process indicator is available upon request.
- ISO 11140 certified inks are water based and non-toxic.
- All materials comply with international standards.
- Opening direction is clearly marked.
- Validated production process.
- Each roll is wrapped in a protective film and placed in the packing carton to ensure protection against dust even after opening of the box.
- Suitable to be used with all types of sealing machines.

SRS - STERILIZATION GUSSETTED REELS 100 M

REF.	SIZE	MIN LOT	PACK
SRS7510	75 mm X 25 mm X 100 m	8 pcs	8 pcs
SRS1010	100 mm X 50 mm X 100 m	4 pcs	4 pcs
SRS1510	150 mm X 65 mm X 100 m	4 pcs	4 pcs
SRS2010	200 mm X 55 mm X 100 m	2 pcs	2 pcs
SRS2510	250 mm X 65 mm X 100 m	2 pcs	2 pcs
SRS3010	300 mm X 80 mm X 100 m	2 pcs	2 pcs
SRS3510	350 mm X 80 mm X 100 m	1 pcs	1 pcs
SRS4010	400 mm X 80 mm X 100 m	1 pcs	1 pcs

STERIDIAMOND GUSSETTED REELS 100 M



“STERIDIAMOND” STERILIZATION FLAT POUCHES

- Sterilization pouches with ‘Chevron’ seal grant a quick and perfect opening.
- Flat pouches imprints are performed outside the packing area, along the welding line, between paper and film, to prevent ink from migrating to the product.
- Transparent green film grants the prompt identification of the content.
- The film is reinforced in order to avoid any tear upon pouch opening.
- High quality “medical grade” paper with 60 or 70 gsm.
- Three-line sealing with strong strength and effective spacing grants a perfect adhesivity during sterilization process and an excellent fibre-free and tear-free peeling during opening.
- Steam and EO process indicators.
- Formaldehyde process indicator is available upon request.
- ISO 11140 certified inks are water based and non-toxic.
- All materials comply with international standards.
- Opening direction is clearly marked.
- Validated production process.
- Flat pouches are packed in plastic bags and placed inside a cardboard box.
- Packs of 100 or 200 pieces are available according to Customer needs.
- Plastic packing ensures protection against dust even after opening of the box.
- Suitable to be used with all types of sealing machines.

SBP - STERIDIAMOND FLAT POUCHES

REF.	SIZE	MIN LOT	PACK	REF.	SIZE	MIN LOT	PACK
SBP7515	75X150 mm.	2000 pcs	2000 pcs	SBP2035	200X350 mm.	1000 pcs	1000 pcs
SBP7520	75X200 mm.	2000 pcs	2000 pcs	SBP2040	200X400 mm.	1000 pcs	1000 pcs
SBP7525	75X250 mm.	2000 pcs	2000 pcs	SBP2045	200X450 mm.	1000 pcs	1000 pcs
SBP7530	75X300 mm.	2000 pcs	2000 pcs	SBP2050	200X500 mm.	1000 pcs	1000 pcs
SBP7535	75X350 mm.	2000 pcs	2000 pcs	SBP2129	210X290 mm.	1000 pcs	1000 pcs
SBP7540	75X400 mm.	2000 pcs	2000 pcs	SBP2135	210X350 mm.	1000 pcs	1000 pcs
SBP7545	75X450 mm.	2000 pcs	2000 pcs	SBP2525	250X250 mm.	1000 pcs	1000 pcs
SBP7550	75X500 mm.	2000 pcs	2000 pcs	SBP2530	250X300 mm.	1000 pcs	1000 pcs
SBP1015	100X150 mm.	2000 pcs	2000 pcs	SBP2535	250X350 mm.	1000 pcs	1000 pcs
SBP1020	100X200 mm.	2000 pcs	2000 pcs	SBP2540	250X400 mm.	1000 pcs	1000 pcs
SBP1025	100X250 mm.	2000 pcs	2000 pcs	SBP2545	250X450 mm.	1000 pcs	1000 pcs
SBP1030	100X300 mm.	2000 pcs	2000 pcs	SBP2550	250X500 mm.	1000 pcs	1000 pcs
SBP1035	100X350 mm.	2000 pcs	2000 pcs	SBP3030	300X300 mm.	1000 pcs	1000 pcs
SBP1040	100X400 mm.	2000 pcs	2000 pcs	SBP3035	300X350 mm.	1000 pcs	1000 pcs
SBP1045	100X450 mm.	2000 pcs	2000 pcs	SBP3040	300X400 mm.	1000 pcs	1000 pcs
SBP1050	100X500 mm.	2000 pcs	2000 pcs	SBP3045	300X450 mm.	1000 pcs	1000 pcs
SBP1055	100X550 mm.	2000 pcs	2000 pcs	SBP3050	300X500 mm.	1000 pcs	1000 pcs
SBP1520	150X200 mm.	1000 pcs	1000 pcs	SBP3055	300X550 mm.	1000 pcs	1000 pcs
SBP1525	150X250 mm.	1000 pcs	1000 pcs	SBP3057	300X570 mm.	1000 pcs	1000 pcs
SBP1530	150X300 mm.	1000 pcs	1000 pcs	SBP3250	320X500 mm.	1000 pcs	1000 pcs
SBP1535	150X350 mm.	1000 pcs	1000 pcs	SBP3658	320X600 mm.	500 pcs	500 pcs
SBP1540	150X400 mm.	1000 pcs	1000 pcs	SBP4040	400X400 mm.	1000 pcs	1000 pcs
SBP1545	150X450 mm.	1000 pcs	1000 pcs	SBP4045	400X450 mm.	1000 pcs	1000 pcs
SBP1550	150X500 mm.	1000 pcs	1000 pcs	SBP4050	400X500 mm.	500 pcs	500 pcs
SBP1555	150X550 mm.	1000 pcs	1000 pcs	SBP4055	400X550 mm.	500 pcs	500 pcs
SBP2020	200X200 mm.	1000 pcs	1000 pcs	SBP4060	400X600 mm.	500 pcs	500 pcs
SBP2025	200X250 mm.	1000 pcs	1000 pcs	SBP5060	500X600 mm.	500 pcs	500 pcs
SBP2030	200X300 mm.	1000 pcs	1000 pcs	SBP6060	600X600 mm.	500 pcs	500 pcs
SBP2035	200X350 mm.	1000 pcs	1000 pcs				

STERIDIAMOND FLAT POUCHES



“STERIDIAMOND” STERILIZATION FLAT POUCHES

- Sterilization pouches with ‘Chevron’ seal grant a quick and perfect opening.
- Flat pouches imprints are performed outside the packing area, along the welding line, between paper and film, to prevent ink from migrating to the product.
- Transparent green film grants the prompt identification of the content.
- The film is reinforced in order to avoid any tear upon pouch opening.
- High quality “medical grade” paper with 60 or 70 gsm.
- Three-line sealing with strong strength and effective spacing grants a perfect adhesivity during sterilization process and an excellent fibre-free and tear-free peeling during opening.
- Steam and EO process indicators.
- Formaldehyde process indicator is available upon request.
- ISO 11140 certified inks are water based and non-toxic.
- All materials comply with international standards.
- Opening direction is clearly marked.
- Validated production process.
- Flat pouches are packed in plastic bags and placed inside a cardboard box.
- Packs of 100 or 200 pieces are available according to Customer needs.
- Plastic packing ensures protection against dust even after opening of the box.
- Suitable to be used with all types of sealing machines.

SBP - STERIDIAMOND FLAT POUCHES

REF.	SIZE	MIN LOT	PACK	REF.	SIZE	MIN LOT	PACK
SBP7515	75X150 mm.	2000 pcs	2000 pcs	SBP2035	200X350 mm.	1000 pcs	1000 pcs
SBP7520	75X200 mm.	2000 pcs	2000 pcs	SBP2040	200X400 mm.	1000 pcs	1000 pcs
SBP7525	75X250 mm.	2000 pcs	2000 pcs	SBP2045	200X450 mm.	1000 pcs	1000 pcs
SBP7530	75X300 mm.	2000 pcs	2000 pcs	SBP2050	200X500 mm.	1000 pcs	1000 pcs
SBP7535	75X350 mm.	2000 pcs	2000 pcs	SBP2129	210X290 mm.	1000 pcs	1000 pcs
SBP7540	75X400 mm.	2000 pcs	2000 pcs	SBP2135	210X350 mm.	1000 pcs	1000 pcs
SBP7545	75X450 mm.	2000 pcs	2000 pcs	SBP2525	250X250 mm.	1000 pcs	1000 pcs
SBP7550	75X500 mm.	2000 pcs	2000 pcs	SBP2530	250X300 mm.	1000 pcs	1000 pcs
SBP1015	100X150 mm.	2000 pcs	2000 pcs	SBP2535	250X350 mm.	1000 pcs	1000 pcs
SBP1020	100X200 mm.	2000 pcs	2000 pcs	SBP2540	250X400 mm.	1000 pcs	1000 pcs
SBP1025	100X250 mm.	2000 pcs	2000 pcs	SBP2545	250X450 mm.	1000 pcs	1000 pcs
SBP1030	100X300 mm.	2000 pcs	2000 pcs	SBP2550	250X500 mm.	1000 pcs	1000 pcs
SBP1035	100X350 mm.	2000 pcs	2000 pcs	SBP3030	300X300 mm.	1000 pcs	1000 pcs
SBP1040	100X400 mm.	2000 pcs	2000 pcs	SBP3035	300X350 mm.	1000 pcs	1000 pcs
SBP1045	100X450 mm.	2000 pcs	2000 pcs	SBP3040	300X400 mm.	1000 pcs	1000 pcs
SBP1050	100X500 mm.	2000 pcs	2000 pcs	SBP3045	300X450 mm.	1000 pcs	1000 pcs
SBP1055	100X550 mm.	2000 pcs	2000 pcs	SBP3050	300X500 mm.	1000 pcs	1000 pcs
SBP1520	150X200 mm.	1000 pcs	1000 pcs	SBP3055	300X550 mm.	1000 pcs	1000 pcs
SBP1525	150X250 mm.	1000 pcs	1000 pcs	SBP3057	300X570 mm.	1000 pcs	1000 pcs
SBP1530	150X300 mm.	1000 pcs	1000 pcs	SBP3250	320X500 mm.	1000 pcs	1000 pcs
SBP1535	150X350 mm.	1000 pcs	1000 pcs	SBP3658	320X600 mm.	500 pcs	500 pcs
SBP1540	150X400 mm.	1000 pcs	1000 pcs	SBP4040	400X400 mm.	1000 pcs	1000 pcs
SBP1545	150X450 mm.	1000 pcs	1000 pcs	SBP4045	400X450 mm.	1000 pcs	1000 pcs
SBP1550	150X500 mm.	1000 pcs	1000 pcs	SBP4050	400X500 mm.	500 pcs	500 pcs
SBP1555	150X550 mm.	1000 pcs	1000 pcs	SBP4055	400X550 mm.	500 pcs	500 pcs
SBP2020	200X200 mm.	1000 pcs	1000 pcs	SBP4060	400X600 mm.	500 pcs	500 pcs
SBP2025	200X250 mm.	1000 pcs	1000 pcs	SBP5060	500X600 mm.	500 pcs	500 pcs
SBP2030	200X300 mm.	1000 pcs	1000 pcs	SBP6060	600X600 mm.	500 pcs	500 pcs
SBP2035	200X350 mm.	1000 pcs	1000 pcs				

STERIDIAMOND FLAT POUCHES



“STERIDIAMOND” STERILIZATION GUSSETTED POUCHES

- Manufactured with a gusset that allows bulky objects to be placed inside
- Sterilization pouches with ‘Chevron’ seal grant a quick and perfect opening.
- Gussetted pouches imprints are performed outside the packing area, along the welding line, between paper and film, to prevent ink from migrating to the product.
- Transparent green film grants the prompt identification of the content.
- The film is reinforced in order to avoid any tear upon pouch opening.
- High quality “medical grade” paper with 60 or 70 gsm.
- Three-line sealing with strong strength and effective spacing grants a perfect adhesivity during sterilization process and an excellent fibre-free and tear-free peeling during opening.
- Steam and EO process indicators.
- Formaldehyde process indicator is available upon request.
- ISO 11140 certified inks are water based and non-toxic.
- All materials comply with international standards.
- Opening direction is clearly marked.
- Validated production process.
- Gussetted pouches are packed in plastic bags and placed inside a cardboard box.
- Packs of 100 or 200 pieces are available according to Customer needs.
- Plastic packing ensures protection against dust even after opening of the box.
- Suitable to be used with all types of sealing machines.

SBS - STERIDIAMOND GUSSETTED POUCHES

REF.	SIZE	MIN LOT	PACK
SBS1030	100X50X300 mm	500 pcs	500 pcs
SBS40	150X50X400 mm	500 pcs	500 pcs
SBS1546	150X50X460 mm	500 pcs	500 pcs
SBS2040	200X55X400 mm	500 pcs	500 pcs
SBS50	200X55X500 mm	500 pcs	500 pcs
SBS2550	250x65x500 mm	500 pcs	500 pcs
SBS55	300x80x500 mm	500 pcs	500 pcs
SBS4258	420X80X580 mm	500 pcs	500 pcs

STERIDIAMOND GUSSETTED POUCHES



“STERIDIAMOND” STERILIZATION SELF-SEALING POUCHES

- Self-sealing pouches do not require the use of thermosealers because they are equipped with a double-sided non-toxic adhesive strip specific for sterilization processes.
- Self-sealing pouches imprints are performed outside the packing area, along the welding line, between paper and film, to prevent ink from migrating to the product.
- Transparent green film grants the prompt identification of the content.
- The film is reinforced in order to avoid any tear upon pouch opening.
- High quality “medical grade” paper with 60 or 70 gsm.
- Three-line sealing with strong strength and effective spacing grants a perfect adhesivity during sterilization process and an excellent fibre-free and tear-free peeling during opening
- Steam and EO process indicators.
- Formaldehyde process indicator is available upon request.
- ISO 11140 certified inks are water based and non-toxic.
- All materials comply with international standards.
- Opening direction is clearly marked.
- Validated production process.
- Standard boxes of 200 pieces including instructions for use.
 Boxes are placed inside a cardboard packaging.
- This kind of packing ensures protection against dust even after opening of the box.

SBA - STERIDIAMOND SELF-SEALING POUCHES

REF.	SIZE	MIN LOT	PACK
SBA60100	60X100 mm	3600 pcs	200 pcs
SBA75250	75X250 mm	1800 pcs	200 pcs
SBA90165	90X165 mm	1800 pcs	200 pcs
SBA09250	90X250 mm	1800 pcs	200 pcs
SBA14260	140X260 mm	1200 pcs	200 pcs
SBA19330	190X330 mm	1200 pcs	200 pcs
SBA19400	190X440 mm	1200 pcs	200 pcs
SBA30370	300X370 mm	800 pcs	200 pcs
SBA30450	300X450 mm	800 pcs	200 pcs
SBA30500	300X500 mm	800 pcs	200 pcs

STERIDIAMOND SELF SEALING POUCHES

