

adhesive wound DRESSINGS

Adhesive wound dressing ULTRA PORE for CATHETER fixation, U shape



Adhesive sterile absorbent U shape dressing, made of soft non-woven material, is designed for fixation of intravenous catheter. Details: Non-woven elastic structure adapts to the contours of the body. The adhesive does not leave any waste. In sterile single package. For single use only. The dressing has an incision in the middle and at the end of the incision, there is a pad that is made from soft viscose, which does not conform to the point of injection. It is simple to use and offers safety and comfort as well as painless removal. When removing the dressing, notice that neither catheters or other devices will not ripped off.

SIZE	ART.	EAN code
6 x 8 cm	157030	4742225009289