

Scalp vein set

Features:

- luer lock
- sterile
- blister
- transparent wing
- shelf life: 5 years

Standards:

- EN ISO 13485
- CE 0123
- class II a

Description of product:


- for single use, non-toxic, pyrogen-free
- venipuncture device intended for the collection of blood, for short-term infusions / transfusions or for the intravenous administration of drugs
- sterilized by ethylene oxide (EO)

Raw material:

- Needle tube: Stainless steel
- Protective cap: Ethylene Polypropylene Copolymer (PP)
- Needle hub with wings: Acrylonitrile Butadiene Styrene (ABS)
- Tubing: Polyvinyl Chloride Polymer
- Male conical fitting: Acrylonitrile Butadiene Styrene (ABS)
- Lubricant: Silicone oil

Packings:

- Blister: 1 piece
- Inner box: 50 pieces
- Outer box: 2000 pieces



Scalp vein set

Standard sizes

Size x length inch	Ext. Dia. x length mm	Tube Length cm	Article-No. [218]
18G x 3/4	1.20 x 19	30	741801D
19G x 3/4	1.10 x 19	30	741900D
20G x 3/4	0.90 x 19	30	742006D
21G x 3/4	0.80 x 19	30	742105D
22G x 3/4	0.70 x 19	30	742204D
23G x 3/4	0.60 x 19	30	742303D
25G x 3/4	0.50 x 19	30	742501D