



Pasūtījumiem un
kontakta informācija:
29178505 ; 20040558



Build Your Imagination

360 patents worldwide; Exported to 140 countries

WOODPECKER



Scan and log on to website
for more information



Find us on:



Find us on:



Find us on:



Guilin Woodpecker Medical Instrument Co., Ltd.

Information Industrial Park, Guilin National High-Tech
Zone, Guilin, Guangxi, 541004 P. R. China
Sales Dept.: +86-773-5873196/2350599
After-sales Service Dept.: 0773-5827898
E-mail: woodpecker4@glwoodpecker.com
Website: <http://www.glwoodpecker.com>

EC REP MedNet EC-Rep GmbH
Borkstrasse 10 · 48163 Muenster · Germany

2024.1.0 Edition

Since **1989**

CATALOGUE

GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.



Tips:

Scaling

G1

Supragingival scaling



G2

Supragingival scaling



G3

Supragingival scaling




G4

Supragingival scaling




G5

Supragingival scaling



G6

Supragingival scaling



G7

Remove dental crown





G8

Remove dental pile, with irrigation



G9

Supragingival scaling



G10

Supragingival scaling



G11

Orthodontic treatment



Periodontics

P1

Periodontal scaling



P3

Periodontal scaling



P4

Periodontal scaling



P3D

Polish root canal surface



P2L

Periodontal scaling



P2R

Periodontal scaling



P2LD

Remove calculus from narrow inter-root spaces (diamond-coated)



P2RD

Remove calculus from narrow inter-root spaces (diamond-coated)



Cavity preparation



SB1

Removing the caries on occlusal surface of teeth and dental neck




SB2

Removing the caries on Mesial surface of the adjacent teeth



SB3

Removing the caries on Distal surface of the adjacent teeth




G31

Preparation of mesiooclusion cavity



G33

Preparation of teeth veneer





SBL

Half ball diamond tip. Set 45° to the left, to remove the carious teeth without damaging the adjacent teeth.



SBR

Half ball diamond tip. Set 45° to the right, to remove the carious teeth without damaging the adjacent teeth.




G30

Preparation between pits and fissures



G32

Preparation of distal and mesial cavity



Shoulder preparation

G20

Abrade dentin



G21

Fine repair and polish the prepared shoulder



Endodontics

E1

120°angle holder, for root canal cleaning





E2

95°angle holder, for root canal cleaning



E3

For the lateral thermal condensation of gutta-percha



E3D

Used to remove the calcification and bad filling material of pulp cavity (diamond-coated)




P4D

Remove calcification in the coronal third of the root canal (diamond-coated)



E4

Remove obstructions and broken instruments in the root canal




E4D

For the retreatment and removal of hard material



E5

Remove root canal fillings and broken instruments in the coronal third





E5D

Remove dental overhangs in the pulp chamber (diamond-coated)



E8

Burs holder, used to expand root canal and grind teeth



E9

Burs holder, used to expand root canal and grind teeth




E10

Used for root canal softly treatment




E10D

Used for root canal retrogression, efficient root apical polishing (diamond-coated)



E11

Used for root canal softly treatment




E11D

Used for root canal retrogression, efficient root apical polishing (diamond-coated)




E14

Remove obstructions and broken instruments in the root canal, with irrigation



E15

Remove root canal fillings and broken instruments in the coronal third, with irrigation



E14D


For the retreatment and removal of hard material, with irrigation



E15D

Remove calcification in root canal, root canal cleaning and enlargement





Tip with excellent tenacity is used for root canal cleaning
File size: 32.5mm±0.5mm
Handpiece size: Φ0.68mm~0.71mm
K-file taper: 2%
Each kit with six pieces,one kit for one model
Contraindications and precautions can be found in the manual.

Scaling (Used for SIRONA)



Supragingival scaling



Supragingival scaling



Supragingival scaling



Supragingival scaling



Supragingival scaling



Supragingival scaling



Remove dental crown



Remove dental pile, with irrigation

Periodontics (Used for SIRONA)



Periodontal scaling



Periodontal scaling



Periodontal scaling



Smoothing restoration overhangs and debride in-depth root (diamond-coated)

Endodontics (Used for SIRONA)



120° angle holder, for root canal cleaning



95° angle holder, for root canal cleaning



For the lateral thermal condensation of gutta-percha



Used to remove the calcification and bad filling material of pulp cavity (diamond-coated)



Remove calcification in the coronal third of the root canal (diamond-coated)



Remove obstructions and broken instruments in the root canal



Remove root canal fillings and broken instruments in the coronal third



Burs holder, used to expand root canal and grind teeth



Used for root canal softly treatment



Used for root canal softly treatment



Remove obstructions and broken instruments in the root canal, with irrigation



Remove root canal fillings and broken instruments in the coronal third, with irrigation

Scaling (Used for KAVO)



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Orthodontic treatment



Supragingival and subgingival scaling



Supragingival and subgingival scaling



Burs holder, used to expand root canal and grind teeth

Periodontics (Used for KAVO)



Subgingival scaling

Used for AMDENT



Supragingival scaling



Supragingival scaling



Supragingival scaling

Bone surgery

US1

Osteotomy

US1L

Osteotomy

US1R

Osteotomy

US2

Osteotomy

US3

Osteotome of precision in the thin regions

US4

Bone remodeling and bone chip harvesting

US5

Bone remodeling and bone chip harvesting

US6

Precise osteotomy

Sinus lifting

UL1

120° Peeling of the sinus membrane

UL2

Peeling of the sinus membrane

UL3

Non-traumatic bone cutting on thin bone region (diamond-coated)

UL4

Micrometric osteotomy near the soft tissue (diamond-coated)*The size of diamond grit is 100µm.

UL5

Peeling of the sinus membrane (angle: 95°)

Endodontic

UE1

Apical root debridement (diamond-coated) *The size of diamond grit is 40µm.

UE2

Soft apical root debridement

UE3

Apical root debridement (diamond-coated) *The size of diamond grit is 40µm.

UE4

Soft apical root debridement

Exelcymosis

UC1

Used for cutting ligaments, minimally damnification exelcymosis

UC3

Used for cutting ligaments, minimally damnification exelcymosis

UC2

Used for cutting ligaments, minimally damnification exelcymosis

Implantation

UI1

Used for implant, with tip ball Φ1.6mm

UI2

Used for implant, with tip ball Φ2mm, irrigation from tip point

UI3

Φ2mm implant site preparation tip,irrigation from tip point.

UI4

Φ3mm implant site preparation tip,irrigation from tip point.

UI5

Φ2mm implant site preparation tip,irrigation from tip point.

UI6

Φ3mm implant site preparation tip,irrigation from tip point.

UI7

Finalize the implant site and bone preparation (diamond-coated)

UI8

For implantation, in plus size, with tip ball Φ1.6mm

UI9

Used for implant, with tip ball Φ2.8mm

Periodontal surgery

UP1

Root scaling

UP2

Root scaling

UP3

Root scaling

UP4

Root surface debridement and root planning (diamond-coated)

*The size of diamond grit is 40µm.

UP5

Perio treatment

UP6

Bone remodeling, root planing in perio operation *The size of diamond grit is 40µm.

UP7

Bone remodeling and root canal enlargement (diamond-coated)

*The size of diamond grit is 40µm.

Ultrasurgery Tip Kit

Basic Kit	Sinus Lifting Kit	Bone Surgery Kit	Exelcymosis Kit	Periodontal Surgery Kit	Endodontic Kit	Implantation Kit
Include: US2/US3 US4/US5/UC1	Include: UL1/UL2 UL3/UL4/UL5	Include: US1/US2 US1L/US1R/US6	Include: UC1/UC2 UC3/UP2/UP3	Include: US4/UP2 UP4/UP5/UP7	Include: UE1/UE2 UE3/UE4/US4	Include: UI1/UI2 UI7/UI8/UI9